



AGENDA

FINAL OF THE POLISH CUP OPEN LONG RANGE SHOOTING

27-28 of SEPTEMBER 2019

(Friday/Saturday)

Prerequisite for participation in the Polish Open Long Range Shooting Cup Final 2019 is participation and obtaining a classified result in one of the earlier rounds of the Polish Open Long Range Shooting Cup 2019.

- **Friday 27th of September 2019**

- training: 9:15 -16:00

- **SATURDAY 28th of September 2019**

- competition: 9:15 - 16:00

- awards ceremony: ~17:30

- **DISTANCE**

Range A - Magnum - distance 600-800-1200m

Range A - Open - distance 300-600-800m

Range B - Standard - distance 300-600-800m

Range B - Historical Rifle - distance 100-200-300m

Range B - Semi -Auto - distance 100-200-300 & 200-300-400m



ADDITIONAL COMPETITION - SUNDAY 29 of SEPTEMBER 2019

The competition is not included in the Cup classification.

Before the rated series at the longest distances in each category, competitors will be given three trial shots.

Players after the trial shots are remain at their shooting positions.

Trial shots will be shown to the players by the judges, based on the photo.

The rated series will begin after the preparation time

- **SUNDAY 29 of SEPTEMBER 2019 R.**

- competition: 9:00 - 15:00
- HAPPY HOURS: 15:00 - 17:00
- awards ceremony: ~17:15
- TELETOMBOLA (prize draw) – ok. 17:30

- **ADDITIONAL DISTANCES**

Range A - Magnum - distance 1200-1400m

Range A - Open - distance 800 – 1000m

Range A - Standard - distance 800 – 1000m

Range B – Historical Rifle - distance 400-600m

Range B - Semi -Auto - distance 400-600m



PREPARATION TIME/NUMBER OF SHOTS

- Magnum **1200m**

(10 shots/10 min)

- Magnum **1400m**

(10 shots/10 min) - 3 trial shots (tape over) time 2 min

- Open/ Standard – **800/1000m**

800m - (10 shots/10 min)

1000m - (10 shots/10 min) - 3 trial shots (tape over) time 2 min

- Historical Rifle – **400/600m**

400m - (10 shots/10 min)

600m - (10 shots/10 min) - 3 trial shots (tape over) time 2 min

- Semi - Auto - **400/600m**

400m - (10 shots/10 min)

600m - (10 shots/10 min) - 3 trial shots (tape over) time 2 min